

Customer:	Location of installation:
_____	_____
_____	_____
_____	_____

Model: _____ SN: _____ Item number: _____
(manual)

The PQ consists of inspections of the correct operation of the cabinet under predefined conditions and procedures. Prerequisites for the PQ are IQ (Installation Qualification) and OQ (Operation Qualification), these must be concluded successfully prior to the initiation of the PQ.

This PQ is intended for the following product series:

BioBlood

Revision: 10/10/2017_001

Person responsible for the cabinet:
Name: _____
Date: _____
Signature: _____
Person responsible for test:
Name: _____
Date: _____
Company: _____
Signature: _____
Person responsible for verification of test:
Name: _____
Date: _____
Company: _____
Signature: _____
Test duration:
Initiation (date/time): _____
Conclusion (date/time): _____

Model: _____ SN: _____

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.
 The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

Measurement - Prerequisites

ID	DESCRIPTION	ACCEPTED	
		YES	NO
P-1	The cabinet must be empty while conducting tests, ie without interior fittings such as drawers, shelves etc. Attachment: Notes:		
P-2	The measurements must be conducted in accordance to IEC 60068-3-5, measured in air with thermocouples or comparable equivalent. Attachment: Notes:		
P-3	The positioning of the sensors in the cabinet must be documented with a sketch and/or a photograph. Attachment: Notes:		

Name:

Signature:

Approved
(Yes / No):

Date:

Conducted by:

Inspected / verified by:

Model:

SN:

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.
 The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

Measurement - Prerequisites

ID	DESCRIPTION	ACCEPTED	
		YES	NO
P-4	Measurements made during the PQ tests must be documented and attached to the PQ. Attachment: Notes:		
P-5	Specify setpoint temperature: _____ °C Specify the ambient temperature: _____ °C Attachment: Notes:		
P-6	Allowed temperature fluctuations - <i>Select the tolerance, according to the model being tested.</i> <i>Find model-specific temperature fluctuations in appendix.</i> Tolerance: +/- _____ K Attachment: Notes:		

Name: _____ Signature: _____ Approved (Yes / No): _____ Date: _____

Conducted by: _____

Inspected / verified by: _____

Model: _____ SN: _____

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.
 The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

Measurement - Temperature stabilization

ID	DESCRIPTION	ACCEPTED	
		YES	NO
P-7	<p>The test is intended to provide substantiation for the temperature stability inside the cabinet during normal operation.</p> <p>The temperature inside the cabinet must be stabilized - where all the points in the working space have reached and maintained the same temperature.</p> <p>When the system is stable, document ordinary operation of the cabinet at the setpoint temperature and ambient temperature specified in P-5.</p> <p>Duration: _____</p> <p>The measurements throughout the operation test, must be documented and attached the PQ.</p> <p>Attachment:</p> <p>Notes:</p>		
P-8	<p>Are the measurements inside the allowed temperature fluctuations specified in P-6 ?</p> <p>Attachment:</p> <p>Notes:</p>		

Name: _____ Signature: _____ Approved (Yes / No): _____ Date: _____

Conducted by: _____

Inspected / verified by: _____

Model: _____ SN: _____

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.
 The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

Measurement - Door opening test

ID	DESCRIPTION	ACCEPTED	
		YES	NO
P-9	<p>The test is intended to provide substantiation for the temperature recovery time inside the cabinet subsequently after a door opening.</p> <p>The temperature inside the cabinet must be stabilized - where all the points in the working space have reached and maintained the same temperature, the setpoint temperature is specified in P-5.</p> <p>When the system is stable, open the door at 90° for 60 seconds.</p> <p>The measurements, throughout the door opening test, must be documented and attached the PQ.</p> <p>Duration: _____</p> <p>Attachment:</p> <p>Notes:</p>		
P-10	<p>Have the setpoint temperature specified in P-5, measured in the absolute centre of the cabinet, been achieved within the set time-frame?</p> <p>Attachment:</p> <p>Notes:</p>		

Name: _____ Signature: _____ Approved (Yes / No): _____ Date: _____

Conducted by: _____

Inspected / verified by: _____

Model: _____ SN: _____

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.
 The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

Measurement - Pull-down

ID	DESCRIPTION	ACCEPTED	
		YES	NO
P-11	<p>The test is intended to provide substantiation for the time it takes for the inside of the cabinet to reach the setpoint temperature specified in P-5. The initial temperature in the working space is the ambient temperature specified in P-5. The temperature inside the cabinet must be stabilized in all points of the working space.</p> <p>When the system is stable. Turn on the power to the cabinet.</p> <p>The measurements, throughout the pull-down test, must be documented and attached the PQ.</p> <p>Attachment:</p> <p>Notes:</p>		
P-12	<p>The time it takes the inside of the cabinet to achieve the setpoint temperature measured in the absolute centre, must not exceed the time-frame specified in the appendix.</p> <p>Duration: _____</p> <p>Have the criteria been met?</p> <p>Attachment:</p> <p>Notes:</p>		

Name:

Signature:

Approved
(Yes / No):

Date:

Conducted by:

Inspected / verified by:

Model:

SN:

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.
 The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

Measurement - Hold-over

ID	DESCRIPTION	ACCEPTED	
		YES	NO
P-13	<p>The test is intended to provide substantiation for the time it takes for the temperature inside the cabinet to reach the terminal temperature specified in the appendix. Ambient temperature and setpoint temperature is specified in P-5.</p> <p>The temperature inside the cabinet must be stabilized - where all the points in the working space have reached and maintained the same temperature throughout, the temperature fluctuations are specified in P-6.</p> <p>When the system is stable, turn off the power to the cabinet.</p> <p>The measurements, throughout the hold-over test, must be documented and attached the PQ.</p> <p>Attachment:</p> <p>Notes:</p>		
P-14	<p>The times it takes the inside of the cabinet to reach the terminal temperature, must at least be the time specified in the appendix.</p> <p>Duration: _____</p> <p>Have the criteria been met?</p> <p>Attachment:</p> <p>Notes:</p>		

	Name: _____	Signature: _____	Approved (Yes / No): _____	Date: _____
Conducted by:	_____	_____	_____	_____
Inspected / verified by:	_____	_____	_____	_____

Model: _____ SN: _____

Deviation Report

Deviations to the criteria of acceptance are to be documented in the deviation report. A separate deviation report shall be made for each deviation. Mark the entry with the relevant "P-ID" specified in the left column in the test specifications.

P-ID: _____

Description of deviation:

Extent to which the deviation has been alleviated:

Additional notes:

Person responsible for test:

Name: _____

Date: _____

Company: _____

Signature: _____

Person responsible for verification of test:

Name: _____

Date: _____

Company: _____

Signature: _____

Model: _____ SN: _____

Approval of test results - Performance Qualification (PQ)

- The steps in the Performance Qualification - PQ were completed with positive results
- The steps in the Performance Qualification - PQ were completed with negative results

ID of steps with negative results: _____

Additional notes:

Person responsible for test

Person responsible for verification of test

Stamp & Signature

Stamp & Signature

Tel.

Tel.

E-mail

E-mail

Location & Date

Location & Date

Model: _____ SN: _____

Appendix:

	Model	Temperature fluctuations	Door opening - recovery time	Pull-down	Hold-over range*	Hold-over
BR	BioBlood					
	500 (Solid door)		3 Minutes	22 Minutes		72 Minutes
	500 (Glass door)		4 Minutes	28 Minutes		42 Minutes
	600D / 600W (Solid door)		3 Minutes	20 Minutes		70 Minutes
	600D / 600W (Glass door)	+/- 2K	4 Minutes	25 Minutes	5°C → 10°C	41 Minutes
	660D / 660W (Solid door)		3 Minutes	20 Minutes		70 Minutes
	660D / 660W (Glass door)		4 Minutes	25 Minutes		41 Minutes
	1270 / 1400 (Solid door)		5 Minutes	23 Minutes		78 Minutes
1270 / 1400 (Glass door)		7 Minutes	29 Minutes		45 Minutes	

* The temperature span between the initial temperature and the terminal temperature in the hold-over test P-13,14

Note:

BR:

Ambient temperature +25°C

Setpoint temperature +5°C

Name:

Signature:

Approved
(Yes / No):

Date:

Conducted by:

Inspected / verified by:

Model: _____

SN: _____



Appendix:

	Model	Temperature fluctuations	Door opening - recovery time	Pull-down	Hold-over range*	Hold-over
BF	BioBlood					
	425		9 Minutes	45 Minutes		55 Minutes
	500		7 Minutes	45 Minutes		55 Minutes
	600D / 600W	+/- 5K	7 Minutes	42 Minutes	-20°C → -10°C	55 Minutes
	660D / 660W		7 Minutes	42 Minutes		55 Minutes
	1270 / 1400		10 Minutes	45 Minutes		58 Minutes
PF	425	+/- 9K	40 Minutes	107 Minutes	-40°C → -10°C	108 Minutes
	600W / 660W	+/- 10K	30 Minutes	57 Minutes	-35°C → -10°C	170 Minutes

* The temperature span between the initial temperature and the terminal temperature in the hold-over test P-13,14

Note:

BF:

Ambient temperature +25°C

Setpoint temperature -20°C

PF (425):

Ambient temperature +25°C

Setpoint temperature -40°C

PF (600W/660W):

Ambient temperature +25°C

Setpoint temperature -35°C

Name:

Signature:

Approved
(Yes / No):

Date:

Conducted by:

Inspected / verified by:

Model: _____

SN: _____

